

MultiTek List of Parts

Part Number	Description
100098	MultiTek, Horizontal, S 115V
100099	MultiTek, Horizontal, N-S-IC 115V
100100	MultiTek, Horizontal, IC 115V
100125	MultiTek, Horizontal, N-S 115V
100126	MultiTek, Horizontal, S 230V
100127	MultiTek, Vertical, N-S 230V
100128	MultiTek, Horizontal, N-S-IC 230V
100129	MultiTek, Vertical, N-S 115V
100130	MultiTek, Vertical, LLS 115V
100131	MultiTek, Vertical, LLS 230V
100132	MultiTek, Horizontal, IC 230V
100135	MultiTek, Horizontal, N 115V
100136	MultiTek, Horizontal, N 230V
100137	MultiTek, Vertical, S 115V
100138	MultiTek, Vertical, S 230V
100139	MultiTek, Horizontal, N-S 230V
100140	MultiTek, Vertical, N 230V
100141	MultiTek, Vertical, N 115V
101407	MultiTek, Horizontal, S-IC 115V
101408	MultiTek, Horizontal, N- IC 115V
101409	MultiTek, Horizontal, S-IC 230V
101410	MultiTek, Horizontal, N- IC 230V
101457	MultiTek, Horizontal, Metrohm IC 115V Bundle
101458	MultiTek, Horizontal, Dionex IC 115V Bundle
101459	MultiTek, Horizontal, Metrohm IC 230V Bundle
101460	MultiTek, Horizontal, Dionex IC 230V Bundle
101461	MultiTek, Horizontal, Metrohm N-S-IC 115V Bundle
101462	MultiTek, Horizontal, Dionex N-S-IC 115V Bundle
101463	MultiTek, Horizontal, Metrohm N-S-IC 230V Bundle
101464	MultiTek, Horizontal, Dionex N-S-IC 230V Bundle
101465	MultiTek, Horizontal, Metrohm S-IC 115V Bundle
101466	MultiTek, Horizontal, Dionex S-IC 115V Bundle
101467	MultiTek, Horizontal, Metrohm S-IC 230V Bundle
101468	MultiTek, Horizontal, Dionex S-IC 230V Bundle
101469	MultiTek, Horizontal, Metrohm N-IC 115V Bundle
101470	MultiTek, Horizontal, Dionex N-IC 115V Bundle
101471	MultiTek, Horizontal, Metrohm N-IC 230V Bundle
101472	MultiTek, Horizontal, Dionex N-IC 230V Bundle
	Options
	Upgrade Kits
101439	MultiTek IC Upgrade Kit
101440	MultiTek Nitrogen Upgrade Kit 115V
101441	MultiTek Nitrogen Upgrade Kit 230V
101442	MultiTek Sulfur Upgrade Kit
101445	MultiTek LLS Upgrade Kit
	Computer

101413	MultiTek PC Bundle
	Router
1700-000-013	115V Ethernet Router 8-port / PAC
1700-000-010	230V Ethernet Router 8-port / PAC
	Sample Introduction
	Liquids
99230-115V	Model 735 Controlled Rate Sample Drive, Horizontal
99230-230V	Model 735 Controlled Rate Sample Drive, Horizontal
99230V-115V	Model 735 Controlled Rate Sample Drive, Vertical
99230V-230V	Model 735 Controlled Rate Sample Drive, Vertical
33127-4-115V	Model 748 Robotic Liquid Autosampler, MiniPal, Vertical
33127-4-230V	Model 748 Robotic Liquid Autosampler, MiniPal, Vertical
33127-6-115V	Model 748 Robotic Liquid Autosampler, MiniPal, Horizontal
33127-6-230V	Model 748 Robotic Liquid Autosampler, MiniPal, Horizontal
102672	Model 758 Robotic Liquid Autosampler, Horizontal, 100V
102671	Model 758 Robotic Liquid Autosampler, Horizontal, 115V
102670	Model 758 Robotic Liquid Autosampler, Horizontal, 230V
102669	Model 758 Robotic Liquid Autosampler, Vertical, 100V
102668	Model 758 Robotic Liquid Autosampler, Vertical, 115V
102667	Model 758 Robotic Liquid Autosampler, Vertical, 230V
	748 Options
33120-10	Syringe Block Kit 10ul - included with purchase of Model 748 (choice of one of 10ul or 25ul)
33120-25	Syringe Block Kit 25ul - included with purchase of Model 748 (choice of one of 10ul or 25ul)
33120-100	Syringe Block Kit 100ul
	Gas/LPG
99087-1-115V	Model 734-1A Automated Gas/Liquefied Gas Sampling Inlet System
99087-1-230V	Model 734-1A Automated Gas/Liquefied Gas Sampling Inlet System
99087-2-115V	Model 734-2A Automated Gas Sampling Inlet System
99087-2-230V	Model 734-2A Automated Gas Sampling Inlet System
99087-3-115V	Model 734-3A Automated Liquefied Gas Sampling Inlet System
99087-3-230V	Model 734-3A Automated Liquefied Gas Sampling Inlet System
99088-1-115V	Model 734-1 Gas/Liquefied Gas Sampling Inlet System
99088-1-	Model 734-1 Gas/Liquefied Gas Sampling Inlet System

230V	
99088-2-115V	Model 734-2 Gas Sampling Inlet System
99088-2-230V	Model 734-2 Gas Sampling Inlet System
99088-3-115V	Model 734-3 Liquefied Gas Sampling Inlet System
99088-3-230V	Model 734-3 Liquefied Gas Sampling Inlet System
	Gas Loop Sample Sizes
21059-0.5CC	Loop Sample S/S 1/8" Treat - 0.5 CC
21059-1CC	Loop Sample S/S 1/8" Treat - 1 CC
21059-2CC	Loop Sample S/S 1/8" Treat - 2 CC
21059-3CC	Loop Sample S/S 1/8" Treat - 3 CC
21059-4CC	Loop Sample S/S 1/8" Treat - 4 CC
21059-5CC	Loop Sample S/S 1/8" Treat - 5 CC
21059-10CC	Loop Sample S/S 1/8" Treat - 10 CC
21059-15CC	Loop Sample S/S 1/8" Treat - 15 CC
21059-20CC	Loop Sample S/S 1/8" Treat - 20 CC
21059-25CC	Loop Sample S/S 1/8" Treat - 25 CC
21059-50CC	Loop Sample S/S 1/8" Treat - 50 CC
21059-70CC	Loop Sample S/S 1/8" Treat - 70 CC
21059-75CC	Loop Sample S/S 1/8" Treat - 75 CC
	LPG Loop Sample Sizes
21023-5uLT	Loop Sample S/S 1/16 Treated - 5 uL
21023-10uLT	Loop Sample S/S 1/16 Treated - 10 uL
21023-20uLT	Loop Sample S/S 1/16 Treated - 20 uL
21023-25uLT	Loop Sample S/S 1/16 Treated - 25 uL
21023-30uLT	Loop Sample S/S 1/16 Treated - 30 uL
21023-40uLT	Loop Sample S/S 1/16 Treated - 40 uL
21023-50uLT	Loop Sample S/S 1/16 Treated - 50 uL
21023-250uLT	Loop Sample S/S 1/16 Treated - 250 uL
	Solids
101719	Model 737-M Temperature Programmable Inlet With Turbo **
	** Requires vacuum pump, P/N 80.24.401 or 80.24.402
80.24.401	Vacuum Pump 115V
80.24.402	Vacuum Pump 230V
99091-0-115V	Model 740 Multi-Matrix Inlet System
99091-0-230V	Model 740 Multi-Matrix Inlet System
99091-2-115V	Model 740 Multi-Matrix Inlet System, for use with Constant Service Model 734

99091-2-230V	Model 740 Multi-Matrix Inlet System, for use with Constant Service Model 734
	Consumables
36117	Syringe 25uLLeap Contip
36118	Syringe, 10ul
37708	Std. - H2O 25PPM Install.
37709	Std. - H2O 5PPM Install.
69090	Septa 12 mm Tef-BAC (Bulk)
70002	Septa, 12.5 mm (12pk)
70047	O-ring, .9891D x .070 S70
70049	O-ring, 1.614D x .070 S70
70061	Pyrotube Liq. Direct Inj
70142	O-ring 1.9251D x .103 S70
70144	O-ring 2.4251D x .103 S70
70150	O-ring, 2.3641D x .070 E70
70173	O-ring, .7991D x .103 S70
71029	Quartz Chips (50Gm / Pk) (2Ln/Pk)
71033	Gas Dryer Drierite / Mol. Sieve
71050	O-ring, .9841D x .139 S70
71051	O-ring .8641D x .070 S70
71069	O-ring, .4261D x .070 S70
71091	O-ring, .426D x .070 Kalrez
71092	O-ring .2391D x .070 Kalrez
71231	O-ring, 1.4871D x .103 N70
37012	Standard Kit NS-Pyr Dbt Iso
37030	Standard Kit NS-Pyr Dms0 H2O
37042	Standard Kit NS-Pyr Dbt Iso (L)
37045	Standard Kit NS-Pyr Dms0 H2O (L)
37043	Standard Kit NS-Pyr Dbt Iso (M)
37046	Standard Kit NS-Pyr Dms0 H2O (M)
37044	Standard Kit NS-Pyr Dbt Iso (H)
37047	Standard Kit NS-Pyr Dms0 H2O (H)
37048	Standard Kit Fluoride Iso
37049	Standard Kit Fluoride H2O
37050	Standard Kit Lls Dbt. Pyr. Iso
37299	Special Order Standard Kit
90073	Syringe, 100ul
101474	Pyrotube Kit: Direct Inject Pyrotube 70061 without Ceramic Insert
101475	Pyrotube Kit: 71323 Multimatrix w/ 734 Support
101476	Pyrotube Kit: 71303 Low Level Sulfur
101478	LEAP Injector Kit 25ul
101479	LEAP Injector Kit 100ul
101496	Kit, O-rings and Filters, Sulfur Reaction Chamber, MultiTek
	758 Consumables
102706	Solvent Vial Caps, 758 Autosampler
102707	Solvent Vial Septa, 758 Autosampler
102708	20mL Vial - SOLVENT VIAL 10pk (100 pieces each)
102709	10mL Vial - WASTE VIAL 10pk (100 pieces each)
102752	Kit, 758, 1 Waste Vial Including Septum + Cap
32057	Vials, Screw Top, Autosampler 100 / Pkg.

36117	Syringe 25uL Leap Contip
90073	Syringe 100ul LEAP
	Spares
16248	Wrist Strap, Anti-Static
17168	Reaction chamber sleeve
18027-1	Clamp. Spring S/S Pyro-Hld
18038	Cover Plate 714
25531	Bkt. Furn. Ceramic Close
30000	Union, Brass Bulkhead 1/8"
30004	Tee, Run, Brass, Female 1/8"
30008	Ferrule Brass 1/8" Set
30018	Union, Br. Bulkhead 1/4"
30029	Ferrule, 1/4" vespel
30060	Tubing Copper, 1/8" od. x .03Wall G
30066	Tubing, teflon 1/8" white
30069	Gas Line Copper 6Ft.
30080	Tubing S/S 1/16 x .040 id.
30081	Tube, 1.4 Teflon, 3/16ID
30099	Thermocouple adaptor fit.
30189	Union Reducer S/S 1/4"-1/16
30426	Fit. 1/8" Plug Brass
30583	Fit. React. Chamb. Inlet
63030	Tubing Shrink 1/8" Black
63031	Tubing Shrink 1/4" Black
70008	Elbow Glass 18/9 Socket x 90
70030-2	Photomultiplier Tube, Nitrogen
70032	Filter. Band Pass 1.25" Dia
70058	Membrane Dryer
70090	PM Tube, R6095HA, Sulfur
70151	TEC, Sealed Assy, .2STG. 12V
70185	Gasket PVC PMT Fl. Ch.
71011	Dry Caps
71078	Ozone cell assby, MD D
71093-1	Lens Collim U. V.
71098	Filter Interference for 90180
71111	In-Line Filter (nitrogen only)
71137-2	Scrubber Cat/Dri 1/4" -1/8" (N/Ns)
79127	Cerafelt (2" x 2")
81045	Boot Fit. Light Leak 1/8"
81050	Boot Light Leak 1/4"
81100	1 Pr. Cera Furn End 1/2" T
82026	Clamp 18 Ball Joint
86021	Dual Thermocouple R6000
86078	Heated Transf Line Assy.
101477	Coking/Sooting Kit: N and NS 70061 Pyrotube
101480	Preventative Maintenance 1 Year LLS or Sulfur Only
101481	Preventative Maintenance 1 Year NS or N Only
101484	Preventative Maintenance 2 Year N Only
101485	Preventative Maintenance Halides 1 Year
101486	Preventative Maintenance Halides 2 Years
101488	740 Service Kit

101489	Nitrogen Moisture Kit
101490	Monthly Maintenance Kit
101491	Coking/Sooting Kit: S 70061 Pyrotube
	758 Spares
102710	Cable, 758, MultiTek/9000 Interface
102711	Adapter, 758, Vertical Injector Port
102712	Screw Kit, 758, Vertical Plate
102713	Screw Kit, 758, Horizontal Plate
102724	Kit, 758, Fuses 230V
102725	Kit, 758, Fuses 115V
102730	Plug, 758, Plunger
102731	Block, 758, Syringe
102733	Tray, 758, 110 Position for 2mL Vials
102747	Kit, 758, Standard Tools
102748	Packaging, 758
102749	Stickers, 758, Antifriction
102751	Kit, 758, 4 Solvent Vials Including Septa + Cap
16350	Rs-232/485 Serial Cable Db9 - Db9
81022	Power Cord 5-15P C13 115V
81025	Power Cord Shuko C13 230V
	Kits
101474	Pyrotube Kit: Direct Inject Pyrotube 70061 without Ceramic Insert
101475	Pyrotube Kit: 71323 Multimatrix w/ 734 Support
101476	Pyrotube Kit: 71303 Low Level Sulfur
101477	Coking/Sooting Kit: N and NS 70061 Pyrotube
101491	Coking/Sooting Kit: S 70061 Pyrotube
101478	LEAP Injector Kit 25ul
101479	LEAP Injector Kit 100ul
101480	Preventative Maintenance 1 Year LLS or Sulfur Only
101481	Preventative Maintenance 1 Year NS or N Only
101484	Preventative Maintenance 2 Year N Only
101485	Preventative Maintenance Halides 1 Year
101486	Preventative Maintenance Halides 2 Years
101488	740 Service Kit
101489	Nitrogen Moisture Kit
101489	Nitrogen Moisture Kit