

	Item No.	Description
<b>Systems</b>		
	15-0500-000	AQUATek 100 Liquid Autosampler, universal voltage with 5mL PEEK sample loop. Including 100 vial (40mL) capacity, 2 position variable dosing standard valves. Includes chiller tray, blank water reservoir, start up kit, required software for Stratum and Velocity operation, and 1 year warranty. If installing to a Velocity XPT with serial number US07320009 or earlier p/n 14-8781-190 D.C. Output board will be required for proper operation. Pricing available on Velocity XPT spare parts list.
<b>Accessories</b>		
	21-0700-000	Recirculating bath, refrigerated, Model 2050W. 120V, 60Hz, 11 amps, compact size - 12-1/2" x 8-3/4" footprint
	21-0701-000	Recirculating bath, refrigerated, Model 2050W. 220V, 60Hz, 11 amps, compact size - 12-1/2" x 8-3/4" footprint
	14-4362-000	Installation kit for 2050W water bath
	15-0528-100	Upgraded Kit, A100 pH Sensor
<b>Sample Loops</b>		
	15-0545-067	Sample Loop, PEEK, 5ml A-100
	15-0545-167	Sample Loop, PEEK, 25ml A-100
	15-0545-267	Sample Loop, PEEK, 20ml A-100
	15-0545-367	Sample Loop, PEEK, 10ml A-100
<b>Electronics</b>		
	15-0551-000	Concentrator Adapter Cable 25 to 9 pin
	14-4383-028	Switch, Power, 10 AMP Filtered use Delta
	14-5665-034	Fuse, Slow Blow, 5AMP, (115V), 5 X 20 MM
	14-8781-490	PCB, DC Control
	14-9593-090	Power Supply, 24 VDC
	14-9639-090	Power Supply 25W
	15-0504-090	PCB, AT-100 Temp Board
	15-0380-190	PCBA, Motor Control, A100
<b>Fittings</b>		
	12-0044-016	Ferrule set, 1/8", brass, set
	12-0069-016	Nut, 1/8", brass
	12-0070-016	Tee, 1/8", Brass
	12-0074-016	Union, bulkhead, 1/8" brass
	12-0075-016	Conn, Male, 1/8" Tube to 1/8" NPT
	14-0053-016	Elbow, male, 1/8", NPT, brass
	14-0056-000	1/8" NPT to 1/8" Compression
	14-0241-016	Ferrule, 1/16", SS
	14-0243-116	Nut, 1/16", Short Gold plated
	14-3123-016	Ferrule, 1/16, Upchurch ETFE (081-042)
	14-7038-016	Bushing, Exit, 1/8" Tefzel, Short
	14-7671-016	Ferrule, Super Flangeless, 1/16", PEEK
	14-7695-016	Nut, PEEK, Natural for 1/16" OD Tubing
	14-8470-016	Plug Nut, 1/4-28, Tefzel, Flat Bottom

	14-8820-016	Nut, 1/16", Short, PEEK
	14-9778-016	Nut 1/8", 1/4 - 28, PEEK, Long
	14-7201-009	Washer, TEFLON
<b>Glassware</b>		
	14-2333-024	Sparger, 25mL fritless, 1/2" (glassware only)
	14-2334-024	Frit Sparger, 25mL, 1/2" neck, Tekmar
	14-2336-024	Fritless Sparger, 5mL, 1/2" Tekmar
	14-2337-024	Frit Sparger, 5mL, 1/2" Neck Tekmar
	14-4006-024	Sparger, 5mL frit w/ side arm (glassware only)
	14-4007-024	Sparger, 25mL frit w/ side arm (glassware only)
	14-4007-224	Sparger kit, 25 mL frit w/ side-arm, Velocity XPT
	14-4826-224	Sparger kit, 25 mL fritless w/ side-arm, Velocity XPT
	14-8921-000	Sparger kit, 5mL fritted
	14-8922-000	Sparger kit, 25 mL fritted
	14-8923-000	Sparger kit, 5 mL fritless
	14-8924-000	Sparger kit, 25 mL fritless
	15-0201-024	Vessel, 15mL Amber IS
	15-0202-079	Bushing, Peek W/ O-Ring, IS Vessel
<b>Heaters</b>		
	14-7806-000	Assy, RTD (24") Temp Sensor for Soil Valve, Mount
	15-0561-020	Assy, Water Reservoir/Heater Universal Voltage
<b>Miscellaneous</b>		
	14-7664-100	Pressurized Reservoir
	14-8632-073	Clamp, Water Heater
	14-9525-300	Rinse Station
	14-9620-294	Drain Tray
	14-9786-079	Vial Stripper Adapter
	15-0457-053	2-Stage Sample Needle
	15-0463-180	Carousel, 100 Position, 40ml Vial
	15-0500-017	Box and Foam
	15-0500-075	A-100 Software & Manual CD
	15-0500-KIT	Kit Box
	14-9266-400	Carousel, Drive Assy, AQUATek 100
	15-0503-000	Pull Handle
	15-0554-073	Clamp - Cable, Adjustable
<b>Tubing</b>		
	14-0583-002	Tubing, 1/8" OD, .062" ID, Natural FEP (per foot)
	14-6934-002	PEEK Tubing, 1/16" OD, Natural (per foot)
	14-7630-002	Tubing, 1/8" OD, Red, Translucent (per inch)
	14-7634-002	Tubing, 1/8" OD, Blue, Translucent (per inch)
	14-8595-002	Tubing, 1/8", Peek, (per foot)
	14-8927-002	PEEK Tubing, 1/16" OD, Green (per foot)
	14-9642-002	Tubing, Wire Reinforced (Clear), ea= 1ft
	15-0272-002	Line Assy, I.S., w/ Minstac Connector
	15-0507-000	Service Loop Tubing Assy
<b>Valves</b>		

	14-3940-000	Pressure gauge, 0-30 psi
	15-0112-100	3 POS Manifold (w/ 1-plug), Internal Standard Valve AQUATek 100
	15-0236-050	Valve, Is Atomx
	15-0271-000	Mounting Kit, IS Valve
	15-0485-050	Valve 24VDC, SMC
	15-0806-200	Regulating, Self Relieving, 0-15psi