

Automated Softening point RB36 5G - 115 / 230V

Part Number	Description
V25001	RB36 5G 115/230V 50/60HZ
	Accessories
K02075	PRINTER 230V 40 COLUMNS
K02076	PRINTER 115V 40 COLUMNS
Q048040	TEST BALL FOR DETECTION CONTROL
Q2500120	POURING PLATE BITUMEN
V02303	PS400 : TEMPERATURE CHECK BOX 0 TO 400°C
V250100	HOT COATING KIT EN1871
V250102	BALL COUPLING CAGES RB36
	Consumable for 1 year
E00904	STIRRING MAGNET D 9 L 45MM
H008A10	DELAYED FUSE 5X20 10A
H008A6.3	DELAYED FUSE 5X20 6A3
H06315	PROTISTOR FUSE 10X38 30A
L001104	RB36 5G PROBE (E2PROM)
M00413	RB36 BEAKER STANDARD
Q2500120	POURING PLATE BITUMEN
V048100	RING & BALL TEST HOLDER
V048101	PACK OF 10 SHOULDERED RINGS
V048105	PACK OF 10 BALLS DIAM 9.52
	Service spare parts for 2 years
C00409	TRIAC 40A/400V
E02014M	FAN 24VCC WITH CONNECTIONS
E02022M	FAN 24VCC 16CFM WITH CONNECTIONS
	All Consumable
E00904	STIRRING MAGNET D 9 L 45MM
H008A10	DELAYED FUSE 5X20 10A
H008A6.3	DELAYED FUSE 5X20 6A3
H06315	PROTISTOR FUSE 10X38 30A
L001104	RB36 5G PROBE (E2PROM)
M00413	RB36 BEAKER STANDARD
M004170	800ML BEAKER EN1871
Q2500120	POURING PLATE BITUMEN
R00232-10	PACK OF 10 BALLS D15 EN1871
U00951	HEAT PROTECTION GLOVES
V048100	RING & BALL TEST HOLDER
V048101	PACK OF 10 SHOULDERED RINGS
V048101-5	SET 5 PACK 10 RINGS V048101
V048105	PACK OF 10 BALLS DIAM 9.52
V048105-5	SET 5 PACK 10 BALLS V048105

V04819	SPECIAL TEST HOLDER EN1871
V250101	RING KIT EN1871
	All Service spare parts
A090164	HEATING ELEMENTS (SHAPE COIL)
C00409	TRIAC 40A/400V
D01127	FLAT KEYPAD 5G SERIES
E015116	TESTED STEPPING MOTOR
E02014M	FAN 24VCC WITH CONNECTIONS
E02022M	FAN 24VCC 16CFM WITH CONNECTIONS
K01068M	GRAPHIC SCREEN BACKLIGHT LED TESTED
K01069	POWER SUPPLY 200W 115/230V 24V
Q048040	TEST BALL FOR DETECTION CONTROL
X21002M	SUBSET SCREEN/KEYB.WITH LABEL
Y04811	INFRA-RED TRANSMITER BOARD EQUIPPED
Y20202	HEATING REGULATION BOARD
Y21001M	CPU BOARD EQUIPPED
Y21003M	LCD JUNCTION BOARD TESTED
Y25001M	TESTED MANAGING BOARD RB36 5G
Y25002	RB36 5G BALL DETECTION BOARD