

<b>MultiTek S-IC</b>	
Part Number	Description
	<b>IC Ready - (No IC included)</b>
101407	MultiTek, Horizontal, S-IC 115V
101409	MultiTek, Horizontal, S-IC 230V
	<b>IC Bundles</b>
101465	MultiTek, Horizontal, Metrohm S-IC 115V Bundle
101467	MultiTek, Horizontal, Metrohm S-IC 230V Bundle
101466	MultiTek, Horizontal, Dionex S-IC 115V Bundle
101468	MultiTek, Horizontal, Dionex S-IC 230V Bundle
	<b>Options</b>
	<b>Computer</b>
101413	MultiTek PC Bundle
	<b>Routers</b>
1700-000-013	115V Ethernet Router 8-port / PAC
1700-000-010	230V Ethernet Router 8-port / PAC
	<b>Upgrade Kits</b>
101440	MultiTek Nitrogen Upgrade Kit 115V
101441	MultiTek Nitrogen Upgrade Kit 230V
	<b>Accessories</b>
	<b>Liquid Sample Accessories</b>
102672	Model 758 Robotic Liquid Autosampler, Horizontal, 100V
102671	Model 758 Robotic Liquid Autosampler, Horizontal, 115V
102670	Model 758 Robotic Liquid Autosampler, Horizontal, 230V
33127-6-115V	Model 748 Robotic Liquid Autosampler, MiniPal, Horizontal
33127-6-230V	Model 748 Robotic Liquid Autosampler, MiniPal, Horizontal
	<b>100 uL syringe recommended for Halides analysis. Either 10 or 25 uL recommended for S analysis (only one block included with purchase of Model 748)</b>
33120-10	Syringe Block Kit 10ul
33120-25	Syringe Block Kit 25ul
33120-100	Syringe Block Kit 100ul
	<b>Gas/LPG Sample Accessories</b>
99087-1-115V	Model 734-1A Automated Gas/Liquefied Gas Sampling Inlet System
99087-1-230V	Model 734-1A Automated Gas/Liquefied Gas Sampling Inlet System
99087-2-115V	Model 734-2A Automated Gas Sampling Inlet System
99087-2-230V	Model 734-2A Automated Gas Sampling Inlet System
99087-3-115V	Model 734-3A Automated Liquefied Gas Sampling Inlet System
99087-3-230V	Model 734-3A Automated Liquefied Gas Sampling Inlet System
99088-1-115V	Model 734-1 Gas/Liquefied Gas Sampling Inlet System
99088-1-230V	Model 734-1 Gas/Liquefied Gas Sampling Inlet System
99088-2-115V	Model 734-2 Gas Sampling Inlet System
99088-2-230V	Model 734-2 Gas Sampling Inlet System
99088-3-115V	Model 734-3 Liquefied Gas Sampling Inlet System
99088-3-230V	Model 734-3 Liquefied Gas Sampling Inlet System

	<b>Gas Loop Sample Sizes</b>
21059-0.5CC	Loop Sample S/S 1/8" Treat - 0.5 CC
21059-1CC	Loop Sample S/S 1/8" Treat - 1 CC
21059-2CC	Loop Sample S/S 1/8" Treat - 2 CC
21059-3CC	Loop Sample S/S 1/8" Treat - 3 CC
21059-4CC	Loop Sample S/S 1/8" Treat - 4 CC
21059-5CC	Loop Sample S/S 1/8" Treat - 5 CC
21059-10CC	Loop Sample S/S 1/8" Treat - 10 CC
21059-15CC	Loop Sample S/S 1/8" Treat - 15 CC
21059-20CC	Loop Sample S/S 1/8" Treat - 20 CC
21059-25CC	Loop Sample S/S 1/8" Treat - 25 CC
21059-50CC	Loop Sample S/S 1/8" Treat - 50 CC
21059-70CC	Loop Sample S/S 1/8" Treat - 70 CC
21059-75CC	Loop Sample S/S 1/8" Treat - 75 CC
	<b>LPG Loop Sample Sizes</b>
21023-5uLT	Loop Sample S/S 1/16 Treated - 5 uL
21023-10uLT	Loop Sample S/S 1/16 Treated - 10 uL
21023-20uLT	Loop Sample S/S 1/16 Treated - 20 uL
21023-25uLT	Loop Sample S/S 1/16 Treated - 25 uL
21023-30uLT	Loop Sample S/S 1/16 Treated - 30 uL
21023-40uLT	Loop Sample S/S 1/16 Treated - 40 uL
21023-50uLT	Loop Sample S/S 1/16 Treated - 50 uL
21023-250uLT	Loop Sample S/S 1/16 Treated - 250 uL
	<b>Solids/Liquids Sample Accessories</b>
99091-0-115V	Model 740 Multi-Matrix Inlet System
99091-0-230V	Model 740 Multi-Matrix Inlet System
99230-115V	Model 735 Controlled Rate Sample Drive, Horizontal
99230-230V	Model 735 Controlled Rate Sample Drive, Horizontal
101719	737-M Temperature Programmable Inlet **
	<b>** Requires vacuum pump, P/N 80.24.401 or 80.24.402</b>
80.24.401	Vacuum Pump 115V
80.24.402	Vacuum Pump 230V
	<b>Recommended Startup Spares and Consumables:</b>
	<b>Standards</b>
37708	Std. - H2O 25PPM Install.
37709	Std. - H2O 5PPM Install.
37012	Standard Kit NS-Pyr Dbt Iso
37030	Standard Kit NS-Pyr DmsO H2O
37042	Standard Kit NS-Pyr Dbt Iso (L)
37045	Standard Kit NS-Pyr DmsO H2O (L)
37043	Standard Kit NS-Pyr Dbt Iso (M)
37046	Standard Kit NS-Pyr DmsO H2O (M)
37044	Standard Kit NS-Pyr Dbt Iso (H)
37047	Standard Kit NS-Pyr DmsO H2O (H)
37048	Standard Kit Fluoride Iso

37049	Standard Kit Fluoride H2O
37299	Special Order Standard Kit
	<b>Other</b>
101475	Pyrotube Kit: 71323 Multimatrix w/ 734 Support
101478	LEAP Injector Kit 25ul
101496	Kit, O-rings and Filters, Sulfur Reaction Chamber, MultiTek
101490	Monthly Maintenance Kit
	<b>Recommended Consumables:</b>
36117	Syringe 25uLLeap Contip
36118	Syringe, 10ul
70002	Septa, 12.5 mm (12pk)
70049	O-ring, 1.614D x .070 S70
70142	O-ring 1.9251D x .103 S70
70144	O-ring 2.4251D x .103 S70
70150	O-ring, 2.3641D x .070 E70
70173	O-ring, .7991D x .103 S70
71029	Quartz Chips (50Gm / Pk) (2Ln/Pk)
71033	Gas Dryer Drierite / Mol. Sieve
71050	O-ring, .9841D x .139 S70
71051	O-ring .8641D x .070 S70
71069	O-ring, .4261D x .070 S70
71091	O-ring, .426D x .070 Kalrez
71092	O-ring .2391D x .070 Kalrez
71231	O-ring, 1.4871D x .103 N70
	<b>Recommended Spares:</b>
101480	Preventative Maintenance 1 Year LLS or Sulfur Only
101485	Preventative Maintenance Halides 1 Year
101486	Preventative Maintenance Halides 2 Years
101488	740 Service Kit
101404	MultiTek IC Cooling Kit