

	Item No.	Description	
Syringe, Valve, and Assemblies			
	14-9784-150	Valve, 7 PORT Spare, (Torch, Phoenix and Apollo,)	
	14-7514-001	Plug, Flat Point Nylon, 1/4-28 X 1/4" LG	
	14-7127-152	Syringe, 25ml with plunger	
	14-8587-30P	Syringe Pumper, Phoenix (NEW)	
	14-9350-001	Screw, Shoulder, 8-32, 3/16 x 3/8, Locking, SS	
	14-8438-179	Syringe Stem	
	14-3661-009	#8x.09 Thick Washer	
	110-252	Knurled Thumb Screw 6-32 x 1/2	
	14-8640-080	Bracket, Syringe Locking	
	14-8661-001	Screw, 6-32 x 1-5 SHCS	
	14-8505-086	Encoder cable assembly, syringe drive pumper	
	14-8455-047	Timing Belt, syringe drive pumper	
	14-8587-30pr	Syringe Pumper, Phoenix (Recertified)	
Glassware and Assemblies in Chemistry Area			
	14-7015-024	Reaction vessel, IC, for Phoenix 8000	
	14-7020-024	Reaction vessel, UV, for Phoenix 8000	
	14-7021-100	UV Reactor Cap Assembly (PEEK)	
	14-7183-000	UV-Lamp with amp connector and O-rings for Phoenix 8000	
	885-379	Cap, sparger TEF	
	14-7022-080	Metal plate for UV reaction vessel cap	
	14-7029-024	Gas/liquid separator for Phoenix 8000	
	14-7033-043	O-ring, sparge tube, in UV reaction vessel	
	050-114	Plug Con 6 CKT Amp 1480704-0	
	885-143	Mist trap for Apollo 9000, DC-180 and Phoenix 8000	
	14-7014-024	Scrubber assembly, chloride, for Apollo 9000 and Phoenix 8000	
	888-448	Tube, scrubber U shape	
	511-738	Pyrex wool 1 gram	
	511-876	Tin 20 mesh 1.5 oz	
	511-895	Copper-20 mesh, 1.5 oz Kit	
	14-8203-000	Kit, Mist Catcher	
	14-7680-024	Fritted sparge tube, IC or UV, for Apollo 9000 or Phoenix 8000	
	100-446	O-Ring, Lamp, UV Reaction Vessel	

	110-236	UV Reactor Screws (teflon cap)	
	14-8809-000	Tool for UV Lamp Lead Pins	
	14-7193-000	Shield, Tinted plastic cover with hinge (Smoked)	
	110-206	UV Reactor Screws (peek cap)	
	14-3823-100	Septum, with 1/8" hole, IC sparger	
Tubing, Fittings, Septa, and Misc in Chemistry Area			
	14-7637-002	Tubing, 1/16OD red, translucent	
	14-7681-002	Tube, TFE 1/16" blue translucent (price per inch)	
	14-0583-002	Tubing, 1/8" OD, .062" ID, Natural FEP (per foot)	
	14-7630-002	Tubing, 1/8" OD, Red, Translucent (per inch)	
	14-7635-002	Tube, 1/8 OD, Green, Translucent	
	14-7634-002	Tubing, 1/8" OD, Blue, Translucent (per inch)	
	14-7633-002	Tubing, 1/8" OD, Yellow, Translucent (per inch)	
	154-424	Tubing, Tygon 1/2" OD (523-359) (price per inch)	
	517-799	Septum, red/white for 1/16th in. tubing	
	517-815	Septa-connector, 1/32" hole, gray, ea (Formerly P/N 517-807, 10 per pack)	
	517-813	Septum, red for 1/8th in. tubing	
	14-0063-043	Septum, plug for mist trap	
	14-8659-016	Nut, Flangless, Short, 1/8, Peek	
	14-8820-016	Nut, 1/16", Short, PEEK	
	14-7035-016	Nut for 1/16" tubing	
	14-7036-016	Bushing, Exit, 1/16" Tezel, Short	
	14-7013-094	Drip tray for Phoenix 8000	
Detector			
	-342322	Cell liner, for Fuji and Apollo IR detectors, 125mm	
	073-022	Cell window, for Fuji and Apollo IR detectors	
	073-023	Cell window o-ring, for Fuji and Apollo IR detectors	
	14-6860-000	Chopper motor/light source assembly for Fuji and Apollo IR detectors	
Plumbing			
	14-0142-016	Connector, Clippard, 10-32-1/8"	
	080-907	Tee 1/8 Polypropylene	
	080-974	Union reducing 1/8 - 1/4 poly	

14-0356-016	Union, bulkhead, 1/8", filter assembly	
090-093	Flow Restrictor 20CC/Min Yellow	
090-097	Flow restrictor, 100mL/min, red dot	
090-099	Flow restrictor, 200CC +10 -5CC 1/8 Union Blu	
14-7026-150	Valve, 2-way air, 24 VDC, for Phoenix 8000	
14-7025-150	Valve, 3-way air, 24 VDC, for Phoenix 8000	
14-7037-016	1/8" Flangless nut	
14-7038-016	Bushing, Exit, 1/8" Tefzel, Short	
090-825	Permeation dryer	
14-7195-000	Flowmeter assembly for Phoenix 8000	
14-4261-016	Cross, 1/8" Union Brass Aps.	

Electronics

14-7009-090	Main CPU board for Apollo 9000 and Phoenix 8000	
14-7008-190	UV power supply for Phoenix 8000	
14-4383-028	Switch, Power, 10 AMP Filtered use Delta	
14-7007-090	Power supply (before S/N 00150001)	
14-7624-090	Power supply (after S/N 00150001)	
14-7717-100	Fan assembly	
14-7682-051	Solid State Relay, 240 VAC, 25A, 0 Volt	
14-7735-000	UV transformer, 115V	
14-7735-100	UV transformer, 230V	

Miscellaneous

14-7045-074	Phoenix 8000 operation manual	
14-8628-000	Filter Set / Maintenance Kit for Whatman 78-40 TOC Gas Generator (40-0014-000)	
14-9515-000	CO2 Filter for Regulator Assy, 14-9512-000 (Replace Annually)	
14-9512-000	Gas Regulator Assy - Includes Gas Regulator and CO2 Filter for carrier gas input into the Apollo.	

DC-183 Boat Section

885-440	PCBA, 183 SSS CNTRL	
---------	---------------------	--