

# BIODIESEL

Part Number	Description
	<b>Biodiesel</b>
CCG4500A	BIODIESEL ALL-IN-ONE SYSTEM ON 7980 120V
CCG4500C	BIODIESEL ALL-IN-ONE SYSTEM ON 7980 230V
	<b>Options:</b>
G4513A	7693A AUTOINJECTOR
G4514A	7693A TRAY, 150 VIAL
G1888A	Headspace Sampler
G1883A	7694E Headspace Sampler
	<b>Accessories:</b>
59.40.200	Oven Exhaust for 6890/7890
	<b>Recommended Startup Spares and Consumables:</b>
	All Spares and Consumables required for a successful startup are included in the startup kit that comes with the instrument
	<b>Recommended 1 Year Consumables:</b>
CCG4500.200	Kit, Consumable Kit Biodiesel All In One
	<b>Recommended 2 Year Spares:</b>
CCG4500.100	Kit, Spare Kit Biodiesel All In One
	<b>Additional Spares:</b>
21040.004	Septa, 9.5 mm (50 pcs)
21040.005	Septa, 11 mm (50 pcs)
21010.194	Seal, TPI (10 pcs)
21032.009	TPI INLET LINER 4MM FOR 0.25/0.32MM COLUMN (10 PCS)
21040.008	FERRULE, GRAPHITE FOR 0.25MM CAP.COL. (10 PCS)
23020.028	Syringe, 5ul ( 2 pcs)
20001.086	Methylheptadecanoate EN 14103
20001.087	C14:0-C24:1 in heptane EN 14103
20001.088	Calibration Set En14105
20001.089	Sample Box Mstfa
20001.090	SAMPLE BOX FORMIC ACID EN14106
20001.091	SAMPLE BOX 2-PROPANOL EN14110
20001.092	SAMPLE BOX METHANOL
20002.047	Sample BioDiesel EN14103
20002.048	SAMPLE MIX BOX EN14105
21052.022	SILICA, DEACTIVATED 3 MTR.
10.73.078	Column, non-Polar 15M, 0.32MM, 0.10UM
21074.032	Column, Innowax 30M, 0.25Mm, 0.25Um Spring-Cage
10.30.201	GC-Base plate for 1 Gasfilter
10.30.210	GC-Triple Gasfilter
27040.100	COLUMN NUT (2 PCS.)