

DHA

Part Number	Description
CCG6510.002A	DHA FAST SYSTEM ON 120V 7890 GC
CCG6510.002C	DHA FAST SYSTEM ON 230V 7890 GC
CCG6510.003A	DHA D 6730 SYSTEM ON 120V 7890 GC
CCG6510.003C	DHA D 6730 SYSTEM ON 230V 7890 GC
CCG6510.006A	DHA D 6729 SYSTEM ON 120V 7890 GC
CCG6510.006C	DHA D 6729 SYSTEM ON 230V 7890 GC
CCG6510.008A	DHA D 6730 COMBI SYSTEM ON 120V 7890 GC
CCG6510.008C	DHA D 6730 COMBI SYSTEM ON 230V 7890 GC
CCG6510.010A	DHA FAST COMBI SYSTEM ON 120V 7890 GC
CCG6510.010C	DHA FAST COMBI SYSTEM ON 230V 7890 GC
CCG6510.013A	DHA D 6729 COMBI SYSTEM ON 120V 7890 GC
CCG6510.013C	DHA D 6729 COMBI SYSTEM ON 230V 7890 GC
CCG6510.014A	DHA D 6733 SYSTEM ON 120V 7890 GC
CCG6510.014C	DHA D 6733 SYSTEM ON 230V 7890 GC
CCG6510.015A	DHA D 6733 COMBI SYSTEM ON 120V 7890 GC
CCG6510.015C	DHA D 6733 COMBI SYSTEM ON 230V 7890 GC
CCG6512.002A	DHA FE SYSTEM ON 120V 7890 GC
CCG6512.002C	DHA FE SYSTEM ON 230V 7890 GC
	6850 GC Systems
ACG6515.002A	DHA FAST SYSTEM ON 120V 6850 GC
ACG6515.002C	DHA FAST SYSTEM ON 230V 6850 GC
ACG6515.003A	DHA D 6730 SYSTEM ON 120V 6850 GC
ACG6515.003C	DHA D 6730 SYSTEM ON 230V 6850 GC
ACG6515.006A	DHA D 6729 SYSTEM ON 120V 6850 GC
ACG6515.006C	DHA D 6729 SYSTEM ON 230V 6850 GC
ACG6515.014A	DHA D 6733 SYSTEM ON 120V 6850 GC
ACG6515.014C	DHA D 6733 SYSTEM ON 230V 6850 GC
	Accessories:
G2880A	6850 Automatic Liquid Sampler
G4513A	7693A AUTOINJECTOR
G4514A	7693A TRAY, 150 VIAL
59.40.200	Oven Exhaust for 6890/7890
	Recommended Startup Spares and Consumables:
	All Spares and Consumables required for a successful startup are included in the startup kit that comes with the instrument
	Recommended 1 Year Consumables:
65186.200	Kit, Consumables DHAfe
65186.400	Kit, Consumables Fast DHA
65186.450	Consumables Fast DHA 6850
65186.800	Kit, Consumables Fast DHA Combi
65186.700	Kit, Consumables DHA D6729
65186.750	Kit, Consumables DHA D6729 (6850)
65186.775	Kit, Consumables DHA D6729 Combi
65186.100	Kit, Consumables DHA D6730
65186.500	Kit, Consumables DHA D6730 Combi
65186.900	Kit, Consumables DHA D6733
65186.975	Kit, Consumables D6733 Combi

DHA

Part Number	Description
	Recommended 2 Year Spares:
65185.850	Kit, Spare Parts (Fast) DHA Combi 7890
65185.800	Kit, Spare Parts DHA-Fe/DHAX-DHA-Combi
	Additional Spares:
10.76.046	ferrule butt connector ID=0,25mm
10.75.046	butt connector (including ferrule) ID=0,25mm
21090.005	Syringe, 0.5uL(2pcs)
20001.114	Sample, n-Paraffin calibration DHAX (5pcs)
20001.115	Sample, Quantitative calibration sample DHAX (5pcs)
20001.116	Sample, n-Paraffin calibration fast DHA (5pcs)
20001.117	Sample, Quantitative calibration sample fast DHA (5pcs)
25115.610	Sample, n-Paraffin calibration DHA Classic (5pcs)
25116.510	Sample, reference RF-A for DHA Classic (5 pcs)
25116.610	Sample, reference RF-C2 for DHAX (5 pcs)
25117.510	Sample, column evaluation DHAX (5pcs)
25115.520	Sample, quantitative calibration sample DHA Classic (5pcs)
20001.220	Sample, Quality Control for fast DHA :Reformer Feed (5 pcs)
20001.221	Sample, Quality Control for fast DHA: Reformate (5 pcs)
20001.222	Sample, Quality Control for fast DHA: FCC Nafta (5 pcs)
20001.223	Sample, Quality Control for fast DHA: Isomerate (5 pcs)
20001.224	Sample, Quality Control for fast DHA: Alkylate (5 pcs)
20001.225	Sample, Quality Control for fast DHA: Gasoline A (MTBE) (5 pcs)
20001.226	Sample, Quality Control for fast DHA: Gasoline B (ETBE+EtOH)(5 pcs)
20001.227	Sample, Quality Control for fast DHA: Gasoline M (EtOH)(5 pcs)
20001.410	Sample, Quality Control DHA D6729: Reformer Feed (5 pcs)
20001.411	Sample, Quality Control DHA D6729: Reformate (5 pcs)
20001.412	Sample, Quality Control DHA D6729: FCC Nafta (5 pcs)
20001.413	Sample, Quality Control DHA D6729: Isomerate (5 pcs)
20001.414	Sample, Quality Control DHA D6729: Alkylate (5 pcs)
20001.415	Sample, Quality Control DHA D6729: Gasoline A (5 pcs)
20001.416	Sample, Quality Control DHA D6729: Gasoline B (5 pcs)
20001.417	Sample, Quality Control DHA D6729: Gasoline M (5 pcs)
20001.420	Sample, Quantitative Calibration DHA D6729 (5 pcs)
20001.421	Sample, n-Paraffin Calibration DHA D6729 (5 pcs)
20001.401	Sample, Quality Control DHA D6730: Reformate (5 pcs)
20001.402	Sample, Quality Control DHA D6730: FCC Naphta (5 pcs)
20001.404	Sample, Quality Control DHA D6730: Alkylate (5 pcs)
20001.405	Sample, Quality Control DHA D6730: Gasoline A (5 pcs)
20001.406	Sample, Quality Control DHA D6730: Gasoline B (5 pcs)
20001.407	Sample, Quality Control DHA D6730: Gasoline M (5 pcs)
20001.409	Sample, Quality Control DHA D6730: Naphta (5 pcs)
20001.121	Sample, D6733 n-Paraffin calibration (5 pcs)
20001.122	Sample, D6733 (512) (5 pcs)
20001.164	Sample, Quality Control for DHA Front End: Reformate (5 pcs)
20001.431	Sample, Quality Control DHA D6733: Reformate (5 pcs)
20001.432	Sample, Quality Control DHA D6733: FCC Naphta (5 pcs)
20001.433	Sample, Quality Control DHA D6733: Isomerate (5 pcs)
20001.434	Sample, Quality Control DHA D6733: Alkylate (5 pcs)
20001.435	Sample, Quality Control DHA D6733: Gasoline A (5 pcs)
20001.436	Sample, Quality Control DHA D6733: Gasoline B (5 pcs)
20001.400	Sample, QC DHA 6730 Reformer feed (5 pcs)
20001.419	Sample, Quality Control DHA D6729 Naphta (5 pcs)
20001.430	Sample, QC DHA D6733 Reformer feed (5 pcs)

DHA

Part Number	Description
21073.038	Column, capillary: Fast DHA L=40m, ID=0.10mm
25190.300	Column, capillary: DHA L=100m, ID=0.25mm
25190.500	Column, capillary DHA D6730 column combination assy L=105m ID 0,25mm
10.70.046	Column, capillary: DB-5 L=15m
21070.037	Column: capillary: DHA L=50m ID=0,25mm
10.74.004	Column, capillary: DHA L=100m, ID=0.25mm on 5" cage
61.21.000	Filter for backflush purge line
21041.001	Ferrule, Graphite for Fast DHA (10 pcs)
10.30.201	GC-Base plate for 1 Gasfilter
10.30.210	GC-Triple Gasfilter
10.32.002	Liner, for DHA-FE Analyzer
21040.004	Septa, 9.5 mm (50 pcs)
21040.005	Septa, 11 mm (50 pcs)
61.30.200	Liner, Liner Assy (pre-frac)
80102.200	Valve, 24V 3-port
27040.100	COLUMN NUT (2 PCS.)