

GAS ANALYZERS

Part Number	Description
Natural Gas Analyzers	
CCG2326A	NGA ISO 6974 SYSTEM ON 120V 7890 GC
CCG2326C	NGA ISO 6974 SYSTEM ON 230V 7890 GC
CCG2327.001A	NGA GPA 2261 1V SYSTEM ON 120V 7890 GC
CCG2327.001C	NGA GPA 2261 1V SYSTEM ON 230V 7890 GC
CCG2327.003A	NGA GPA 2261 3V SYSTEM ON 120V 7890 GC
CCG2327.003C	NGA GPA 2261 3V SYSTEM ON 230V 7890 GC
CCG2328.001A	NGA GPA 2286 EXT. SYSTEM ON 120V 7890 GC
CCG2328.001C	NGA GPA 2286 EXT. SYSTEM ON 230V 7890 GC
6850 GC NGA Systems	
ACG2500A	NGA GPA 2261 SYSTEM ON 120V 6850 GC
ACG2500C	NGA GPA 2261 SYSTEM ON 230V 6850 GC
LPG, Etylene, Propylene Analyzers	
CCG6840A	LPG & PROPYLENE SYSTEM ON 120V 7890 GC
CCG6840C	LPG & PROPYLENE SYSTEM ON 230V 7890 GC
6850 GC LPG, Etylene, Propylene Analyzers	
ACG2410A	LPG & PROPYLENE SYSTEM ON 120V 6850 GC
ACG2410C	LPG & PROPYLENE SYSTEM ON 230V 6850 GC
ACG2400A	S IN ETH.&PROP. SYSTEM ON 120V 6850 GC
ACG2400C	S IN ETH.&PROP. SYSTEM ON 230V 6850 GC
Refinery Gas Analyzers	
CCG2330A	RGA HI-SPEED SYSTEM ON 120V 7890 GC
CCG2330C	RGA HI-SPEED SYSTEM ON 230V 7890 GC
CCG2335A	RGA HI-SPEED+LSV SYSTEM ON 120V 7890 GC
CCG2335C	RGA HI-SPEED+LSV SYSTEM ON 230V 7890 GC
CCG2340A	RGA FAST 1058 SYSTEM ON 120V 7890 GC
CCG2340C	RGA FAST 1058 SYSTEM ON 230V 7890 GC
CCG2345A	RGA FAST 1058+LSV SYSTEM ON 120V 7890 GC
CCG2345C	RGA FAST 1058+LSV SYSTEM ON 230V 7890 GC
Options:	
68100.000	SAMPLING DEVICE WITHOUT VAPORIZER
68100.001	SAMPLING DEVICE INCL.VAPORIZER 110 VOLT
68100.002	SAMPLING DEVICE INCL.VAPORIZER 220 VOLT
68100.003	SAMPLING DEVICE INCL.VAPORIZER 120 VOLT HIGH TEMP
68100.004	SAMPLING DEVICE INCL.VAPORIZER 230 VOLT HIGH TEMP
Accessories:	
Recommended Startup Spares and Consumables:	
All Spares and Consumables required for a successful startup are included in the startup kit that comes with the instrument	
Recommended 1 Year Consumables:	
68486.055	KIT, CONSUMABLES NGA ISO 6974 ON 7890 GC
68486.038	KIT, CONSUMABLES NGA-GPA 2261 GAS/LIQ.
68486.054	KIT, CONSUMABLES NGA GPA 2286
ACG2400.200	KIT, CONSUMABLES S IN ETH.&PROP ANALYZER ON 6850 GC

GAS ANALYZERS

Part Number	Description
CCG2330.200	KIT, CONSUMABLES HI-SPEED RGA
CCG2340.200	KIT, CONSUMABLES FAST RGA ON 7890 GC
	Recommended 2 Year Spares:
68485.055	KIT, SPARE PARTS NGA ISO 6974 ANALYZER
68485.056	KIT, SPARE PARTS NGA GPA 2261 1V
68485.057	KIT, SPARE PARTS NGA GPA 2261 (3VLV)
68485.054	KIT, SPARE PARTS NGA GPA 2286
68485.051	KIT, SPARE PARTS LPG/PROPYLENE ANALYZER
ACG2400.100	KIT, SPARE PARTS S IN ETH.&PROP ANALYZER ON 6850 GC
CCG2330.100	Kit, Spare Parts Hi-Speed Rga
CCG2335.100	KIT, SPARE PARTS HI-SPEED RGA LSV ON 7890 GC
CCG2340.100	KIT, SPARE PARTS FAST RGA ON 7890 GC
CCG2345.100	KIT, SPARE PARTS FAST RGA LSV ANALYZER ON 7890 GC
	Additional Spares:
10.10.161	Inlet Filter, 0.5 micron filter
15.20.007	Press fitting
G1530A.508	Sample loop, 50 ul
G1530A.509	Sample loop, 100 ul
G1530A.502	Gas Sampling Valve Loop 0.25 ml
10.12.006	Rotameter, 0-70ml/min 65mm/nvlv
G1530A.974	Sample Loop 50 uL Sulfinert
G1530A.975	Sample Loop 100uL Sulfinert
G1530A.976	Sample Loop 250ul Sulfinert
99.11.010	Regulator for calibration gas for RGA (29910.020)
99.10.020	Calibration gas for RGA, without regulator (1ltr)
10.73.049	Column, capillary: Pre-column Hi speed L=3m, ID=0.32mm
21071.501	Column, micropacked: HaysepQ 0.25m
21071.502	Column, micropacked: Haysep 1m
21071.503	Column, micropacked: HayesepQ/Molsieve 5A 3m
21071.504	Column, micropacked: Molsieve 13X 2m
84.50.510	Column, packed: Combination Column RGA silcosteel
84.51.000	Column, packed: Pre-column
84.56.000	Column, packed: Molsieve 13X; RGA
84.58.100	Column, packed: Hayesep Q
84.13.000	Column, packed: Molsieve 5A assembly
84.55.000	Column, packed: Hayesep Q; RGA
84.51.010	Column, packed: CW1540 on Por. C Silcosteel PR
84.58.110	Column, packed: Hayesep Q. L=1.2m Silcosteel
18710-20119	FID jet packed
21025.086	Ferrule, for Capillary Column rotary valve Hi-Speed (5 pcs)
10.30.202	GC-Base plate for 2 Gasfilters
10.30.210	GC-Triple Gasfilter
30.30.400	Actuator, repair set, 2 pistons
10.25.005	Actuator, Valco for 10-port Valves
10.25.011	Valve, Air Actuator 6-port (valco)
28431.003	Actuator, for LSV,6,10 port Valves in leftside compartment
10.25.029	Valve, rotary 6-port nitronic 60
10.11.018	TPI solenoid cooling valve
10.11.106	SOLENOID, 5/2 24V FOR INLET BLOCK
10.25.031	VALVE, ROTARY 10-PORT NITRONIC 60
29910.020	TUNING GAS FOR RGA, INCLUDING REGULATOR (1LTR)
27040.100	COLUMN NUT (2 PCS.)