

SIMDIS

Part Number	Description
CCG2123.894A	SIMDIS D 2887 SYSTEM ON 120V 7890 GC, includes AC FastSimdis and AC8634
CCG2123.894C	SIMDIS D 2887 SYSTEM ON 230V 7890 GC, includes AC FastSimdis and AC8634
CCG2123.884A	SIMDIS D 3710 SYSTEM ON 120V 7890 GC
CCG2123.884C	SIMDIS D 3710 SYSTEM ON 230V 7890 GC
CCG2123.886A	SIMDIS HT 750 SYSTEM ON 120V 7890 GC
CCG2123.886C	SIMDIS HT 750 SYSTEM ON 230V 7890 GC
CCG2123.890A	SIMDIS D 5442 WAX SYSTEM ON 120V 7890 GC
CCG2123.890C	SIMDIS D 5442 WAX SYSTEM ON 230V 7890 GC
CCG2123.891A	SIMDIS D 7096 SYSTEM ON 120V 7890 GC
CCG2123.891C	SIMDIS D 7096 SYSTEM ON 230V 7890 GC
CCG2123.892A	SIMDIS D 7213 SYSTEM ON 120V 7890 GC
CCG2123.892C	SIMDIS D 7213 SYSTEM ON 230V 7890 GC
CCG2120AA	AC 8612 SYSTEM ON 120V 7890 GC
CCG2120AC	AC 8612 SYSTEM ON 230V 7890 GC
	Dual Channel Systems
CCG2125.894A	DC SIMDIS D 2887 SYSTEM ON 120V 7890 GC
CCG2125.894C	DC SIMDIS D 2887 SYSTEM ON 230V 7890 GC
CCG2125.884A	DC SIMDIS D 3710 SYSTEM ON 120V 7890 GC
CCG2125.884C	DC SIMDIS D 3710 SYSTEM ON 230V 7890 GC
CCG2125.886A	DC SIMDIS HT 750 SYSTEM ON 120V 7890 GC
CCG2125.886C	DC SIMDIS HT 750 SYSTEM ON 230V 7890 GC
CCG2125.890A	DC SIMDIS D 5442 WAX SYSTEM ON 120V 7890 GC
CCG2125.890C	DC SIMDIS D 5442 WAX SYSTEM ON 230V 7890 GC
CCG2125.891A	DC SIMDIS D 7096 SYSTEM ON 120V 7890 GC
CCG2125.891C	DC SIMDIS D 7096 SYSTEM ON 230V 7890 GC
	Additional Channels (incl. inlet and detector)
	● Allows for additional channel only when first channel is also AC Simdis channel.
CCG2127A.894	ADDITIONAL SIMDIS D 2887 CHANNEL ON 7890 GC, includes AC FastSimdis and AC8634
CCG2127A.884	ADDITIONAL SIMDIS D 3710 CHANNEL ON 7890 GC
CCG2127A.886	ADDITIONAL SIMDIS HT 750 CHANNEL ON 7890 GC
CCG2127A.890	ADDITIONAL SIMDIS D 5442 WAX CHANNEL ON 7890 GC
CCG2127A.891	ADDITIONAL SIMDIS D 7096 CHANNEL ON 7890 GC
CCG2127A.892	ADDITIONAL SIMDIS D 7213 CHANNEL ON 120V 7890 GC
	● Channels can only run simultaneously if both channels are the same.
	● Column lifetimes might be impacted if maximum temperature is exceeded.
	● Dual channels systems to be configured with additional Auto injector
	Crude Oil Systems
CCG4200A	AC CRUDE OIL ANALYZER SYSTEM ON 120V, INCL DHA FE
CCG4200C	AC CRUDE OIL ANALYZER SYSTEM ON 230V, INCL DHA FE
CCG4200A.001	AC CRUDE OIL SYS. ON 120V, INCL. FAST DHA COMBI
CCG4200C.001	AC CRUDE OIL SYS. ON 230V, INCL. FAST DHA COMBI
CCG4200A.002	AC CRUDE OIL SYS. ON 120V, INCL. DHA D6729 COMBI
CCG4200C.002	AC CRUDE OIL SYS. ON 230V, INCL. DHA D6729 COMBI

SIMDIS

Part Number	Description
CCG4200A.003	AC CRUDE OIL SYS. ON 120V, INCL. DHA D6730 COMBI
CCG4200C.003	AC CRUDE OIL SYS. ON 230V, INCL. DHA D6730 COMBI
	6850 GC Systems
ACG2120AA	AC 8612 SYSTEM ON 120V 6850 GC
ACG2120AC	AC 8612 SYSTEM ON 230V 6850 GC
	Application Kits add an additional mode to an existing channel
	They do not contain a second inlet, detector, or Software
CCG2124A.881	SIMDIS D 2887 APPLIC. ON 7890 GC
CCG2124A.884	SIMDIS D 3710 APPLIC. ON 7890 GC
CCG2124A.886	SIMDIS HT 750 APPLIC. ON 7890 GC
CCG2124A.889	SIMDIS FAST APPLIC. ON 7890 GC
CCG2124A.890	SIMDIS D 5442 WAX APPLIC. ON 7890 GC
CCG2124A.891	SIMDIS D 7096 APPLIC. ON 7890 GC
CCG2124A.892	SIMDIS D 7213 APPLIC. ON 7890 GC
CCG2124A.894	SIMDIS D2887 APPLIC. MODE KIT ON 7890 GC, INCL. AC FAST SIM
	Dual Channel Systems
CCG2126A.881	SIMDIS D 2887 APPLIC. ON 7890 GC DC
CCG2126A.884	SIMDIS D 3710 APPLIC. ON 7890 GC DC
CCG2126A.886	SIMDIS HT 750 APPLIC. ON 7890 GC DC
CCG2126A.889	SIMDIS FAST APPLIC. ON 7890 GC DC
CCG2126A.894	SIMDIS D2887 APPLIC. ON 7890 GC DC, INCL. AC FAST SIMDIS A
	Accessories:
G4513A	7693A AUTOINJECTOR
G4514A	7693A TRAY, 150 VIAL
G2880A	6850 Automatic Liquid Sampler
	Recommended Startup Spares and Consumables:
	All Spares and Consumables required for a successful startup are included in the startup kit that comes with the instrument
	Recommended 1 Year Consumables:
65986.203	Kit, Consumables D2887 Cap (Tpi) On 7890 GC
65986.102	Kit, Consumables D3710 (Tpi) On 7890 GC
65986.350	Kit, Consumables Ht Simdis (Tpi) On 7890 GC
65986.460	Kit, Consumables D7096 Simdis On 7890 GC
65986.465	Kit, Consumables Simdis D7213 On 7890 GC
65986.475	CONSUMABLE KIT D2887 ACCELERATED 7890
65986.455	Kit, Consumables D7096/6850
65986.410	Kit, Consumables AC 8634 /6850
65986.420	Kit, Consumables AC 8612 /6850
	Recommended 2 Year Spares:
65985.010	Kit, Spare Parts Simdis TPI On 7890 GC
	Additional Spares:
23020.028	Syringe, 5ul (2 pcs)
21040.004	Septa, 9.5 mm (50 pcs)

SIMDIS

Part Number	Description
40.01.066	Column cutter,Stainless Steel, High Temp SIMDIS
25901.300	Liner holder PTV
25901.350	Liner, glass for PTV (10pcs)
21010.194	Seal, TPI (10 pcs)
21032.006	TPI liner/seal replacement set (10 pcs)
21032.004	Liner, TPI Inlet 4mm, for 530 um columns (10 pcs)
25650.100	Sample, reference oil, High Temperature (5pcs)
25650.150	Sample, reference oil, ASTM D2887/Fast SIMDIS (5pcs)
25650.250	Sample, reference oil, ASTM D3710 (5pcs)
25950.110	Sample, calibration mix heavy, High Temperature (5pcs)
25950.120	Sample, calibration mix light, High Temperature (5pcs)
25950.150	Sample, calibration mix, ASTM D2887/Fast SIMDIS (5pcs)
25950.200	Sample, calibration mix, extended ASTM D2887 (5pcs)
25950.250	Sample, calibration mix, ASTM D3710 (5pcs)
20001.310	Sample, reference oil, D2887 extended (5pcs)
20001.309	Sample, Graphimetric blend, SimDis HT (5 pcs)
20001.130	Sample, Quality Control SIMDIS 2887 / 8634: Naphtha
20001.131	Sample, Quality Control SIMDIS D2887/ 8634: Jet fuel
20001.550	Sample, Qualitative for Crude Oil Analyzer (5 pcs)
20001.172	Sample, Ref. SD 5442 WAX
20001.553	Sample, Calibration sample, ASTM D7096 (5pcs)
20001.554	Sample, Calibration sample with MTBE, ASTM D7096 (5 pcs)
20001.555	Sample, Calibration light Simdis D7213 (5pcs)
20001.551	Sample, Calibration heavy Simdis D7213 (5pcs)
20001.557	Sample, Reference sample Simdis D7213 (5pcs)
24001.065	Column, capillary, SS: High Temp (2pcs)
24001.070	Column, capillary: SIMDIS D2887 (2 pcs)
24001.075	Column, capillary: extended 2887/AC 8634 (2 pcs)
24001.059	Column, capillary: SIMDIS HT750 Aluclad (2 pcs)
21074.008	Column, capillary: Simdis D2887 on 5" cage (2 pcs)
21073.064	Column, capillary: Extended Fast Simdis D2887 8m (2 pcs)
21073.076	Column, capillary: Simdis D7213 (2 pcs)
21040.009	Ferrule, graphite 1/8" for liner holder PTV SIMDIS (10 pcs)
21040.007	Ferrule, graphite 1/16" x 0.8mm for widebore column (10 pcs)
23020.041	Jet, FID Simdis Capillary, 43mm
24001.006	Jet, FID Simdis Packed/Capillary, 64mm
10.30.201	GC-Base plate for 1 Gasfilter
10.30.210	GC-Triple Gasfilter
27040.100	COLUMN NUT (2 PCS.)