

Automated TAG Flash Point FP56 5G2	
Part Number	Description
V20801	FP56 5G2 115/230V 50/60HZ GAS
	Options and Accessories
K02075	PRINTER 230V 40 COLUMNS
K02076	PRINTER 115V 40 COLUMNS
L00434	ELECTRIC GLOW PLUG FP 5G2
V02301	PS100 : TEMPERATURE CHECK BOX -100 TO 100°C
V074103	TAG CUP
V208101	FP56 5G2 LID
X20406	AIRFLOW PROTECTION
	Consumable for 1 year
L00194	FP 5G2 GLASS TEMPERATURE PROBE
L00419	FP THERM. FLASH POINT DETECTOR
L00433	GAS IGNITOR FP 5G2
L00434	ELECTRIC GLOW PLUG FP 5G2
M00609	SILICON TUBE D3.5X1,5
01000-256-51	CRM FLASH PT TAG KERO 255ML
	All Consumable
H008A3.15	DELAYED FUSE 5X20 3A15
H008A6.3	DELAYED FUSE 5X20 6A3
H06314	THERMAL FUSE SF70
L00194	FP 5G2 GLASS TEMPERATURE PROBE
L00194/3C	FP93 GLASS PROBE 0,40,150°C
L00419	FP THERM. FLASH POINT DETECTOR
L00433	GAS IGNITOR FP 5G2
L00434	ELECTRIC GLOW PLUG FP 5G2
M00609	SILICON TUBE D3.5X1,5
M00610	SILICON TUBE D5.5X3
01000-256-51	CRM FLASH PT TAG KERO 255ML
	Service spare parts for 2 years
C00402	TRIAC 15A/400V TXAL 2215B
D01011	THERMAL PROTECTION
E02022M	FAN 24VCC 16CFM WITH CONNECTIONS
Q2040207	GAS TEST FLAME
Q2060803	PROBE INSULATING GUIDE
	All Service spare parts
C00402	TRIAC 15A/400V TXAL 2215B
D01011	THERMAL PROTECTION
D01125	FP 5G FLAT KEYBOARD
E015116	TESTED STEPPING MOTOR
E02022M	FAN 24VCC 16CFM WITH CONNECTIONS
H06316	PROTISTOR FUSE 10X38 8A
K01052	LCD MONOCHROME SCREEN 1/4 VGA
K01069	POWER SUPPLY 200W 115/230V 24V
L00429	BATH THERMOCOUPLE
M00170	CRYOGENIC SOLENOID VALVE 24V DC
M00203	T-VALVE CHROME PLATED BRASS OLI D3
M00627	GAS HOSE D6X12

Automated TAG Flash Point FP56 5G2	
Part Number	Description
M00807	ROPE CLAMP GAS TUBE
M09806	GAS SOLENOID VALVE FP 5G2
N12002	PILOT FLAME TUBE HOLD SCREW
Q2040107	REST LID
Q2040207	GAS TEST FLAME
Q2040802	SHUTTER TAG
Q2060307	INSULATION SPACER
Q2060803	PROBE INSULATING GUIDE
X20604M	COMPLETE LCD SCREEN
Y20004	2 STEPPING MOTOR BOARD
Y20401	END OF STROKE BOARD
Y20601M	FP ACQUISITION BOARD TESTED
Y20602	LCD JUNCTION BOARD
Y20801	HEATING REGULATION BOARD W/O FUSE
Y21001M	CPU BOARD EQUIPPED