

	<b>Automated Cleveland Flash &amp; Fire Point FP92 5G2</b>
<b>Part Number</b>	<b>Description</b>
V20701	FP92 5G2 115/230V 50/60HZ
V20702	FP92 5G2 ANALYSER WITH SC6
V20703	FP92 5G2 ANALYSER FOR ASPHALTS
	<b>Options and Accessories</b>
K02075	PRINTER 230V 40 COLUMNS
K02076	PRINTER 115V 40 COLUMNS
L00192	FP92 5G2 METAL TEMPERATURE PROBE
Q0420117	CLEVELAND CUP WITHOUT HANDLE
V02303	PS400 : TEMPERATURE CHECK BOX 0 TO 400°C
V030103	OVERFLOW CUP WITH HANDLE CLEVELAND
V20002	GAS PRESSURE REGULATION 3,7 KPA MAX.
V20710	SC6 COC SAMPLE CHANGER
V20711	COC ASPHALT TESTING KIT
	<b>Consumable</b>
	<b>For standard apparatus:</b>
H008A10	DELAYED FUSE 5X20 10A
H008A6.3	DELAYED FUSE 5X20 6A3
H06315	PROTISTOR FUSE 10X38 30A
H06319	THERMAL FUSE SF152
L00190	FP92 5G2 GLASS TEMPERATURE PROBE
L00192	FP92 5G2 METAL TEMPERATURE PROBE
L00432	GAS IGNITOR CLEVELAND
M00605	GAS HOSE IN. D 1.5 LENGTH 1M
M00760	INSULATION PLATE / CLEV.1P
M00770	CLV INSULATION PLATE
Q0420117	CLEVELAND CUP WITHOUT HANDLE
Q0420118	OVERFLOW CUP CLEVELAND WITHOUT HANDLE
Q2070912	REMOVABLE RAKE / ASPHALT OPTION
V030101	CLEVELAND TEST CUP WITH HANDLE STANDARD
V030103	OVERFLOW CUP WITH HANDLE CLEVELAND
V207105	CLEVELAND COMPLETE METAL PLATE
01000-292-51	CRM FLASH PT COC OIL 255ML
	<b>Consumable</b>
	<b>For sample changer:</b>
M00768	SC6 INSULATION PLATE
Q0420117	CLEVELAND CUP WITHOUT HANDLE
Q0420118	OVERFLOW CUP CLEVELAND WITHOUT HANDLE
	<b>Consumable</b>
	<b>For asphalt testing option:</b>
Q2070912	REMOVABLE RAKE / ASPHALT OPTION
	<b>Service spare parts for 2 years</b>
A090164	HEATING ELEMENTS (SHAPE COIL)
C00409	TRIAC 40A/400V
D01125	FP 5G FLAT KEYBOARD
D01126	SC6 FRONT PANEL / KEYBOARD
E015116	TESTED STEPPING MOTOR
E015121	MOTOR KIT FOR SC6 SERIAL NUMBER ->200
E02022M	FAN 24VCC 16CFM WITH CONNECTIONS

	<b>Automated Cleveland Flash &amp; Fire Point FP92 5G2</b>
<b>Part Number</b>	<b>Description</b>
H06315	PROTISTOR FUSE 10X38 30A
J01080	LINK CABLE RS 9 PINS / 3 METERS
J01096	DETEC.FLAME TEST WIRE CLEVELAND
K01052	LCD MONOCHROME SCREEN 1/4 VGA
K01053	CONVERTER+CORD FOR K01052
K01069	POWER SUPPLY 200W 115/230V 24V
L00432	GAS IGNITOR CLEVELAND
M00164	CONNECTOR FOR SV COIL T22
M00203	T-VALVE CHROME PLATED BRASS OLI D3
M00627	GAS HOSE D6X12
M00765	HARD SCREEN PROTECTION
M00770	CLV INSULATION PLATE
M00807	ROPE CLAMP GAS TUBE
M09805	GAS SOLENOID VALVE CLV
Q0420117	CLEVELAND CUP WITHOUT HANDLE
Q0420118	OVERFLOW CUP CLEVELAND WITHOUT HANDLE
Q2020701	PLEXI PROTECTION
R00419	TEETH BELT MXL 110 TEETH
R00421	SYNCHRONUS BELT 175 TEETH
V207110	HEATING BLOCK
X20604M	COMPLETE LCD SCREEN
X20701M	CLV MOTOR ASSEMBLY TESTED
X20706	IONIZ.RING STANDARD CLV
X20708	IONIZ.RING ASPHALT V207
Y20004	2 STEPPING MOTOR BOARD
Y20007	FIRE DETECTION BOARD
Y20201	END OF STROKE BOARD
Y20202	HEATING REGULATION BOARD
Y20301	SAMPLE CHANGER CONTROL PANEL
Y20601M	FP ACQUISITION BOARD TESTED
Y20602	LCD JUNCTION BOARD
Y21001M	CPU BOARD EQUIPPED