

Automatic Cold Filter Plugging Point FPP Analyzer	
Part Number	Description
V22101	FPP 5GS 100-240V 50/60HZ
V22101ML	FPP 5GS 100-240V 50/60HZ ML
	Options and Accessories
K02075	PRINTER 230V 40 COLUMNS
K02076	PRINTER 115V 40 COLUMNS
V02301	PS100 : TEMPERATURE CHECK BOX -100 TO 100°C
V21001	VACUUM CONTROL KIT 300MM
V220100	12V DC POWER SUPPLY OPTION
V221103	MANUAL CPP KIT FOR FPP5G/5GS
	Consumables for 1 year
L00186	FPP 5G TEMPERATURE PROBE
M004232	400ML BEAKER
M004234	FPP 5G TEST TUBE
M004235	FPP 5G PIPETTE
M00591	MEMBRANE FOR MICRO-PUMP M00583
M00607	Silicone Tube D7.5x5mm Lg1m
M01028	PACK OF 10 GASKETS 5.28 CREP CFPP
M01029	PACK OF 10 GASKETS 12.42 CREP CFPP
M01056	SET OF 10 VITON GASKET 5.7X 1.9
Q0330303	PACK OF 10 FILTER GAUZE
Q2210401	FPP INSULATING RING AND SPACERS
V211103	FPP 5G FILTRATION ASSEMBLY
V033102	COMPLETE FILTER UNIT
V211105	TUBE CONNECTOR WITH JOINTS
01000-309-51	CRM CFPP TRANSITION BLEND 255ML
01000-309-52	CRM CFPP WINTER BLEND 255ML
	Service spare parts for 2 years

L00189	PT100 SENSOR D3 BRASS CGRID
M00169	3 WAYS SOLENOID VALVE 24V DC
M00583	MICRO-PUMP 12VCC
	All Consumables
L00186	FPP 5G TEMPERATURE PROBE
L00503	ASTM 6C THERMOMETER
L00504	ASTM 5C THERMOMETER
M004232	400ML BEAKER
M004234	FPP 5G TEST TUBE
M004235	FPP 5G PIPETTE
M00591	MEMBRANE FOR MICRO-PUMP M00583
M00592	VALVE FOR MICRO-PUMP M00583
M00607	Silicone Tube D7.5x5mm Lg1m
M01011	SEALING RING D5,70X1,9
M01028	Pack of 10 gaskets 5,28 inlet filter CFPP
M01029	PACK OF 10 GASKETS 12.42 CREP CFPP
M01056	SET OF 10 VITON GASKET 5.7X 1.9
M01104	CORK RING D 43,8
M01105	CORK DISK D 45 X 6
M01126	CONICAL CORK D34-D29 H34 DRILL6,7
Q0330303	PACK OF 10 FILTER GAUZE
Q0330304	TEST FILTER HEAD MAIN PART
Q0330307	KNURLED UPPER NUT
Q0330308	CLAMP FOR CONICAL RING
Q0330309	SPHERICAL RING
Q0330311	METAL GAUZE HOLDER
Q0330312	FILTER FASTENING NUT
Q2110117	BEAKER STOPPER
Q2210401	FPP INSULATING RING AND SPACERS
V033102	COMPLETE FILTER UNIT
V211103	FPP 5G FILTRATION ASSEMBLY
V211105	TUBE CONNECTOR WITH JOINTS
01000-309-51	CRM CFPP TRANSITION BLEND 255ML
01000-309-52	CRM CFPP WINTER BLEND 255ML

	All Service spare parts	
D01013M	THERMAL PROTECT 60°C WIRED	
D01127	FLAT KEYPAD 5G SERIES	
E02022M	FAN 24VCC 16CFM WITH CONNECTIONS	
H00429	SHIPMENT LOCKING DEVICE	
J02038	POWER CONNECTOR + SWITCH	
K01068M	GRAPHIC SCREEN BACKLIGHT LED TESTED	
K01069	POWER SUPPLY 200W 115/230V 24V	
L00189	PT100 SENSOR D3 BRASS CGRID	
L00754	OPTICAL DETECTOR FPP 5G	
M00169	3 WAYS SOLENOID VALVE 24V DC	
M00230	DISCHARGE LIMITER D2 G1/8	
M003191	EXHAUST SILENCER M 1/8	
M004233	FLASK D56 X H100	
M00583	MICRO-PUMP 12VCC	
M00591	MEMBRANE FOR MICRO-PUMP M00583	
M00592	VALVE FOR MICRO-PUMP M00583	
M005987	PUMP SILENT	
M01011	SEALING RING D5,70X1,9	
M01022	VITON JOINT DI 29,30XDT 3,60	
V221104	FPP 5GS COMPLETE JACKET	
X21002M	SUBSET SCREEN/KEYB.WITH LABEL	
Y21001M	CPU BOARD EQUIPPED	
Y21002M	COOLING CONTROL BOARD TESTED	
Y21003M	LCD JUNCTION BOARD TESTED	
Y21101M	FPP5G ACQUISITION BOARD TESTED	