

GUM TESTER

Part Number	Description
09151-000-01	Gum test apparatus 3-post-unit, 230V, 50/60 Hz. air jet evaporation method
09152-000-01	Gum test apparatus 3-post-unit, 115V, 50/60 Hz. air jet evaporation method
09171-000-01	Gum test apparatus, 3-post unit, 230V 50/60 Hz air/steam jet evaporation method, with built-in steam superheater
09172-000-01	Gum test apparatus, 3-post unit, 115V 50/60 Hz air/steam jet evaporation method, with built-in steam superheater
	Accessories:
09181-000-00	Steam generator mod. BM 2090 - 9kw, 400V / 50/60Hz. Switchable in 3 steps (3+3+3 kw), overpressure max. 9 bar, vapour consumption max.13 kg/h, boiler housing made of stainless steel, fittings made of brass. Bleeding through ball cock R1/2"
09182-000-00	Feedwater vessel BM030 made of stainless steel,content approx. 30 ltr., incl. high pressure pump for warm water operation max. 70°C
09191-000-00	Compressor for air evaporation test / Gumwill be supplied complete with water trap module and pressure reducer Sucking capacity: 200 l/min. Eff. supply volume at 3.5 bar: 140 l/min. Max. pressure: 7 bar Tank volume: 20 l - stationary
343-039	Protective tube, ni. / Gumtest
502-047	Thermometer 3C ASTM E1, -5+400 C/ Gum
	Spare parts:
304-027	Pressure gauge 0-1,0 bar / Gumtest
343-034	Conical air/steam evaporation tube
412-080	Pressure reducer w.o. manometer
501-012	Test jar for Gum-Test
502-047	Thermometer 3C ASTM E1, -5+400 C/ Gum
613-050	Solid state relay 25A
617-053	Main switch/green (mains supply) 230V
617-159	Overtemperature protection 280°C HGT
617-160	Temperature regulator 2-fold/DMP48AS/HGT917

GUM TESTER

Part Number	Description
620-060	Heating tape 230V 1100W/Gums
620-061	Heating plate 230V 420W/Gum, AP-Semi
620-112	Heating plate 230V 520W Gumtest
628-007	Fuse 10 A - MT
628-028	Fuse 16 A - T
628-033	Fuse holder FEB 031. 1401
636-114	Pt100 to 300°C/HGT917