

HDA 620 DISTILLATION (manual)	
Part Number	Description
06201-000-00	HERZOG Manual Distillation Apparatus HDA 620 STANDARD Version, 230 V, 50/60 Hz
06201-000-01	HERZOG Manual Distillation Apparatus HDA 620 STANDARD Version, 115 V, 50/60 Hz.
06201-000-06	HERZOG Manual Distillation Apparatus HDA 620 with BLOWER and FIRE EXTINGUISHER 230 V 50/60 Hz.
06201-000-07	HERZOG Manual Distillation Apparatus HDA 620 with BLOWER and FIRE EXTINGUISHER 115 V 50/60 Hz.
	Accessories:
110-165	Bath control unit 230V/50/60Hz for HDA 620
110-166	Bath control unit 115V/50/60Hz for HDA 620
	Spare parts:
101-075	Ball cock, brass nick. R3/8" Reid Semi/HDA620
101-129	Ball valve G 1/4" brass ni.-plated
108-213	Push button switch 617-101 with mounting tool 617-105
108-289	Cover for condenser bath HDA 620
121-009	O-ring 15,60 x 1,78 Viton
244-280	Adjustment gauge f. flask heater helix and holder for cleaning set
301-016	Distillation drip tray for measuring cylinder
305-050	Ceran board with 25 mm central opening
305-051	Ceran board 32 mm central opening
305-052	Ceran board with 38 mm central opening
305-053	Ceran board with 50 mm central opening
412-212	Hard wood pad
413-040	Water gauge
413-043	Stopper, Silicone, 15mm OD with hole 6.5 (flask to condenser)
413-125	Centering stopper, Teflon, f. 125 ml distillation flask 510-011 w. NS19/26 ground neck, acc.to ASTM D 86/DIN 51751
413-127	Centering stopper, Teflon, for 200 ml flask 510-012 with NS24/29 acc. to ASTM D 1078
413-220	Centering stopper, Teflon, for 125 ml distillation flask 510-005 straight neck flask acc. to ASTM D 86 Fig.4
413-225	Centering stopper, Teflon, for 200 ml flask 510-006 with streight neck acc. to ASTM D 86
502-020	Thermometer 7C ASTM E1, -2+300°C/Dist.manual
502-021	Thermometer 8C ASTM E1, -2+400°C/Dist.manual
502-022	Thermometer 7F ASTM E1, +30+580°F/Dist.manual
502-023	Thermometer 8F ASTM E1, +30+760°F / Dist.manual
501-160	Drop glass D20x180 mm for vapor thermometer parking
510-005	Distillation flask 125 ml with straight neck and with thermometer position mark
510-006	Distillation flask 200 ml
510-009	Distillation flask 100 ml, straight neck
510-011	Distillation flask 125 ml with mark and NS 19/26 ground neck and with thermometer position mark

HDA 620 DISTILLATION (manua	
Part Number	Description
510-012	Distillation flask, 200 ml with ground neck NS24/29 and with thermometer position mark
511-001	Measuring cylinder 100 ml / Dist.man.
607-218	Heating Control Board HDA 620
613-065	Solid-State-Relay 50 A for distillation heater
617-053	Main switch/green (mains supply) 230V
625-111	S-Trafo for HDA620 115/208/230V-24V/1200VA
628-006	Fuse 0.5 A - MT
628-007	Fuse 10 A - MT
628-018	Fuse 0.05 A - MT
628-028	Fuse 16 A - T
639-008	Fan for oven / HDA 620
645-046	Heating Helix for flask heater, 24V/1000W HDA