

HFP 380 PENSKY MARTENS (ma	
Part Number	Description
03801-000-09	HFP380 Manual Pensky Martens Flash Point Gas Ignition 230 V 50/60 Hz
03801-000-10	HFP380 Manual Pensky Martens Flash Point Gas Ignition, Cooling Blower 230 V 50/60 Hz
03801-000-13	HFP380 Manual Pensky Martens Flash Point Gas Ignition 115 V 50/60 Hz
03801-000-14	HFP380 Manual Pensky Martens Flash Point Gas Ignition, Cooling Blower 115 V 50/60 Hz
	Spare parts:
103-071	Test cover PM Semi gas ignition HFP 380
301-008	Rubber hose/ignitors /flash p.test
302-026	Cup, brass
403-110	Gas connection HFP-Semi/manual - Gas ignition
502-003	Thermometer 9C ASTM E1, -5+110 C / PM
502-004	Thermometer 10C ASTM E1, +90+370 C / PM
502-099	Thermometer 88C ASTM E1, +10+200 C /PM
502-005	Thermometer 9 F ASTM E1, +20+230 F /PM
502-006	Thermometer 10F ASTM E1,+200+700 F /PM
615-061	LED red
617-003	Micro switch
617-036	Rocker switch + retaining s pring / COC Semi
617-053	Main switch/green (mains supply) 230V
618-137	Motor with gear, 24V-3,9W/280RPM for HFP380
618-138	Motor with gear, 24V-3,9W/54RPM for HFP380
627-098	Cooling fan for MP 329 + HDA 626/627/628
628-023	Fuse 6.3 A - MT
628-007	Fuse 10 A - MT
628-008	Fuse 0.1 A - MT
628-018	Fuse 0.05 A - MT
645-031	Bottom heater 230V 450W HFP36x/38x/329
645-032	Wall heater 500W MP329/330, HFP36x
645-059	Wall heater 115 V / 500W / HFP 380
645-071	Heating plate 115V/400W HFP380