

HRB 754 RING & BALL	
Part Number	Description
07541-000-01	Automatic Ring & Ball Tester HRB 754 230V 50/60Hz
07541-000-02	Automatic Ring & Ball Tester HRB 754 115V 50/60Hz
07549-110-70	Optional blower module for HRB 754 For extended bath temperature to +250°C (To be ordered with the instruments)
	Accessories:
101-214	Citizen Printer DotMatrix 100-230V
101-212	Ink ribbon black/red Citizen
101-185	Paper roll for Matrix dot printer 101-175/6 + 101-214
626-053	Connection cable PC to printer (CENTRONICS)
101-137	Ink jet printer with German and English manual, f. 100-230 VAC
253-124	Ball coupling cages for test insert HRB 754
253-133	Casting plate for test ring HRB 754/RB36 5G acc. to EN 1427
413-037	Beaker grip tong R+B
636-110	Metal resistance thermometer Pt 100 f. HRB 754, sample
1000-004-001	Pt100 probe simulator box 0 to +400°C with 6 fixed temp.points, accuracy better 0,2%, with cable
626-217	CAN bus connection cable
619-119	CAN bus terminal load resistor
	Recommended spare parts for appr. 2 years operation:
121-010	O-ring 9,2x2,62 Vito n/MC 452+MC 470
303-042	R&K stirring magnet, TEFLON coated
305-011	Test ring shouldered /R+B-teser acc. To ASTM D 36, E 28, IP 58, ISO 4625, DIN 52011, NFT 66008, EN 1427
305-013	Test ball for R+B tester
326-233	Protection foil for touch panel
504-003	Beaker 800 ml, low form, for R+B and Wilhelmi
501-214	Test tube 10x100mm for HRB 754
628-007	Fuse 10 A - MT
628-020	Fuse 1 A MT
628-023	Fuse 6.3 A - MT
628-025	Fuse 4 A - MT
613-050	Solid state relay 25A
615-104	Voltage regulator 24V TO3
615-108	Laser-Module for HRB 754
616-072	Voltage regulator 78 H 05 SG
617-151	Mains switch 2-pol /230V/50Hz HFP 36x
618-142	Stepper motor for HRB 754
622-012	Rectifier VJ 148 Var o
636-105	Resistance thermometer, glass PT 100 for HRB 754
645-056	Heating 1000 W / HRB 754
	Other spare parts:
121-010	O-ring 9,2x2,62 Vito n/MC 452+MC 470
126-037	Touch Panel with Software / HRB754
110-185	Color touch panel upgrade kit for display 126-035 HRB 754 with serial no. lower than #0048

HRB 754 RING & BALL

Part Number	Description
305-091	Test insert for ring- and ball HRB 754
501-214	Test tube 10x100mm for HRB 754
607-057	Stepper motor/control P.C.B. for stirrer drive
607-097	Thermocouple amplifier P.C.B.
650-143	CPU with software for HRB 754
607-219	Power+Amplifierboard for HRB754
607-224	Detector Board for HRB 754 ball detection
607-226	Interface P.C.B. HRB754
615-104	Voltage regulator 24V TO3
616-072	Voltage regulator 78 H 05 SG
617-151	Mains switch 2-pol /230V/50Hz HFP 36x
618-142	Stepper motor for HRB 754
619-013	Varistor metallic oxide S20 K300
622-012	Rectifier VJ 148 Var o
636-106	Thermocouple/metal/HRB 754
650-143	CPU with software for HRB 754
650-551	Detection board for PT100 programmed