

HVP 906 REID SEMI AUTOMATIC

Part Number	Description
09063-000-00	Digital Reid Semi Automatic Tester complete with 3 test chambers, with 3 pressure sensors and digital displays EN 12, ISO 3007, ASTM D 323 / 4953 DIN 51 754, IP 69 230 V +/- 10%, 50/60 Hz
09064-000-00	Digital Reid Semi Automatic Tester complete with 3 test chambers, with 3 pressure sensors and digital displays EN 12, ISO 3007, ASTM D 323 / 4953 DIN 51 754, IP 69 115 V +/- 10%, 50/60 Hz
	Accessories:
09081-000-00	Air chamber/Reid vapour pressure
09091-000-00	Liquid chamber with 1 opening
502-030	Thermometer -38+50 C:0.1 air chamber / Reid
340-050	Liquid chamber filling device
340-051	Table wrench to screw the bombs
340-052	Connection hose, cpl., REID Semi to testing device for spring gauges
	Recommended spare parts for appr. 2 years operation:
121-012	O-ring 2,2/1mm/for plug Reid S.
121-038	O-ring 2,9/1,78 mm/ Reid S./MC 452
301-026	Seals for bombs/Reid Teflon
340-040	Coupling cpl. with helix/Reid Therm.
340-045	Helix cpl./Reid Thermostat
502-030	Thermometer -38+50 C:0.1 air chamber / Reid
502-034	Thermometer 23C IP +35+42C/Reid
613-002	Relay E 1020 / F 110 24 V
613-019	Control relay 15
613-038	Reed-relay ERID-A-1W 4,1 E
613-050	Solid state relay 25A
617-002	Micro switch
620-049	Heater 230V 1000W /Reid-Semi
628-006	Fuse 0.5 A - MT
628-011	Fuse 6.3 A - F
628-023	Fuse 6.3 A - MT
636-026	Resistance thermometer/Reid S.+MC 625-bath
	Other spare parts:
611-065	Power pack+motor control/Reid Semi
645-039	Heater 115V 1000W /Reid-Semi
618-148	Motor 230V, 50/60Hz, 32W, 2400RPM
618-027	Motor with gear 230/ 240V 50 Hz
611-066	LED-control / Reid Semi
628-007	Fuse 10 A - MT
628-014	Fuse 10.00 A - F
613-022	Switch relay 115V,V23008-A0015-A101
638-003	Float switch cpl./Reid Semi