

	Item No.	Description
Concentrators		
	14-9800-100	StratUm PTC Sample Concentrator 110V. Includes 60" Siltek heated transferline, Siltek tubing and Siltek treated fittings throughout the sample pathway, heated sample mount, water removal system, U-Trap technology, glassware bypass injection port, 5 ml frit sample sparger, Vocarb 3000 trap, #9 trap, spare fuses, instruction manual, and RS-232 interface capability. Includes VOC Teklink software. Requires for operation a PC with CD-ROM drive running Microsoft Windows XP Only. Also requires the appropriate GC or GC/MS interface cable.
	14-9800-1e0	StratUm PTC Sample Concentrator 110V with Guardian Foam Sensor / Foam Eliminator. Includes 60" Siltek heated transferline, Siltek tubing and Siltek treated fittings throughout the sample pathway, heated sample mount, water removal system, U-Trap technology, glassware bypass injection port, 5 ml frit sample sparger, Vocarb 3000 trap, #9 trap, spare fuses, instruction manual, and RS-232 interface capability. Includes VOC Teklink software. Requires for operation a PC with CDROM drive running Microsoft Windows XP Only. Also requires the appropriate GC or GC/MS interface cable. If operating the StratUm PTC without an autosampler a rinse reservoir or rinse tap kit is required.
	14-9800-200	StratUm PTC Sample Concentrator 230V. Includes 60" Siltek heated transferline, Siltek tubing and Siltek treated fittings throughout the sample pathway, heated sample mount, water removal system, U-Trap technology, glassware bypass injection port, 5 ml frit sample sparger, Vocarb 3000 trap, #9 trap, spare fuses, instruction manual, and RS-232 interface capability. Includes VOC Teklink software. Requires for operation a PC with CD-ROM drive running Microsoft Windows XP only. Also requires the appropriate GC or GC/MS interface cable.
	14-9800-2e0	StratUm PTC Sample Concentrator 230V with Guardian Foam Sensor / Foam Eliminator. Includes 60" Siltek heated transferline, Siltek tubing and Siltek treated fittings throughout the sample pathway, heated sample mount, water removal system, U-Trap technology, glassware bypass injection port, 5 ml frit sample sparger, Vocarb 3000 trap, #9 trap, spare fuses, instruction manual, and RS-232 interface capability. Includes VOC Teklink software. Requires for operation a PC with CDROM drive running Microsoft Windows XP only. Also requires the appropriate GC or GC/MS interface cable. If operating the StratUm PTC without an autosampler a rinse reservoir or rinse tap kit is required
	14-9800-075	VOC Teklink Control Software capable of controlling the StratUm PTC with AQUATek 70 / SOLATek 72 autosampler systems
PTC Options		
	14-8417-200	Duet High Throughput Interface. Interfaces two Teledyne Tekmar PTCs with any GC or GC/MS system.
	14-8834-000	Guardian Foam Sensor option
	14-8834-300	Foam Sensor/Eliminator, Stratum
	14-8920-200	Sample Heater Assembly for all glassware types, 110V
	14-8920-300	Assy, Heater, Sample 5/25mL 230VAC

	14-5092-200	Installation Kit, Includes 5ml syringe with luer, 10 ul syringe and sample mount wrench.
	14-7645-000	SOLATek 72/AQUATEk 70 rinse water reservoir kit. Includes 5 liter reservoir, activated carbon filter, copper tubing and fittings.
	14-7644-000	SOLATek 72/AQUATEk 70 external tap rinse kit. Includes activated carbon filter, fittings, copper tubing, saddle valve and pressure regulator assembly.
	14-5745-090	Serial port expansion board for PC (required if your computer has no available serial ports)

Purge & Trap Accessories

	14-8993-000	Cryofocusing Module for StratUm PTC, 110V. Cryofocusing trap for use with capillary columns from 0.53 to 0.20mm ID (non-split applications). Instruction manual and interface to StratUm PTC included. Requires StratUm PTC for operation. If connecting to an Agilent 6890 you must order a 6890 bracket kit P/N 14-7341-280. (For backpressure regulated systems you must purchase an external regulator assembly P/N 14-3938-000)
	14-8993-100	Cryofocusing Module for StratUm PTC, 230V. Cryofocusing trap for use with capillary columns from 0.53mm to 0.20mm ID (non-split applications). Instruction manual and interface to StratUm PTC included. Requires StratUm PTC for operation. If connecting to an Agilent 6890 you must order a 6890 bracket kit P/N 14-7341-280. (For backpressure regulated systems you must purchase an external regulator assembly P/N 14-3938-000)
	14-7341-280	HP 6890 GC mounting bracket
	14-3938-000	External regulator assembly. Required when interfacing to any type of Mass Selective Detector and a cryofocusing module. (Must bypass flow controller when using MSD) when flow rates are under 1ml/min.
	14-4814-000	LN(2) tank adapter kit
	14-8856-000	Assy, 10mL Antifoam 1520

PTC Glassware Options

	14-8921-000	Sparger kit, 5mL fritted
	14-8922-000	Sparger kit, 25 mL fritted
	14-8923-000	Sparger kit, 5 mL fritless
	14-8924-000	Sparger kit, 25 mL fritless
	14-8925-100	kit, 5mL needle sparge, Stratum
	14-8926-100	kit, 25mL needle sparge, Stratum

Cables

	14-1197-039	Cord, Power, 115V 15 A, European Recept
	14-6689-086	Interface, HP 6850, 6890, 7890 GC
	14-2991-000	Interface, HP 5890 GC
	14-4652-086	Interface, HP 5890 w/5970 MSD and Unix-B or MS-DOS software, HP 5890/5971/5972 MSD and Unix-B or MS-DOS software and HP 5890/5989 MS Engine
	14-4188-086	Interface, HP 5890 w/5970 MSD and Unix or Pascal based software. Note: This cable requires the HP A111 (HP P/N 05990-60111) or A211 (HP P/N 05990-60211) accessory card and internal accessory cable (HP P/N 05987-60158) if using the Pascal software. If using Unix software the internal accessory cable (HP P/N 05987-60158) only is necessary.
	14-2993-000	Interface, HP 5995/96/85/87/88/92 GC/MS with HP-1000/RTE GC/MS Software.

14-7396-086	Interface, Varian 3800/3900
14-3233-000	Interface, Perkin-Elmer 8000 Series/Autosystem, Clarus.
14-5397-086	Interface, two Tekmars on one Perkin-Elmer 8000 Series/Autosystem
14-4610-086	Interface, Shimadzu GC 14A/15A/17A GC 14A w/QP 1000 EX MSD and GC 14A w/QP 2000 MSD, also includes 1/16" to 2mm unions, needs Shimadzu P/N 221-34280-91
14-4938-086	Interface, Carlo Erba Mega and Vega Series/Fisons 8000
14-3147-000	Interface, general purpose/HNU 301/321/421 (Valve driver option necessary from HNU.)
14-7397-086	Trace 2000/Trace GC Ultra Interface cable w/ 1/16" to 2mm unions
14-7397-186	Trace 2000/Trace GC Ultra Interface cable
14-3233-200	PE Clarus 500 Interface Cable
14-9754-086	Shimadzu 2010-2014 Interface Cable
14-8477-286	Cable, Archon-Tekmar Interface
14-5028-039	230V Power Cable